

Product Datasheet

Human UbcH5d protein (orb49714)

Description	Recombinant Human Ubiquitin-Conjugating Enzyme E2 D4 is produced by our E.coli expression system and...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-GST
Note	For research use only
Application notes	GST tag at the N-terminus.
Uniprot ID	Q9Y2X8
Expiration Date	6 months from date of receipt.