

Product Datasheet

Human UBE2J2 protein (orb49717)

Description	Recombinant Human Ubiquitin-Conjugating Enzyme E2 J2 is produced by our E.coli expression system and...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.2.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-GST
Note	For research use only
Application notes	GST tag at the N-terminus.
Uniprot ID	Q8N2K1
Expiration Date	6 months from date of receipt.