

Product Datasheet

Human ATG3 protein (orb49722)

Description	Recombinant Human Ubiquitin-like-conjugating Enzyme ATG3/ATG3 is produced by our E.coli expression s...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Application notes	T7 tag at the N-terminus, 6His tag at the C-terminus.
Uniprot ID	Q9NT62
Expiration Date	6 months from date of receipt.