

## Product Datasheet

### Mouse FGF acidic protein (orb49774)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Recombinant Mouse/Rat Fibroblast growth factor 1/Fibroblast Growth Factor Acidic is produced by our ...   |
| <b>Endotoxins</b>      | Less than 0.1 ng/ug (1 EU/ug)   |
| <b>Form/Appearance</b> | Lyophilized from a 0.2 um filtered solution of PBS,PH7.4.   |
| <b>Storage</b>         | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| <b>Tag</b>             | tag free  |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>P61148</b>   |
| <b>Expiration Date</b> | 6 months from date of receipt.  |