


## Product Datasheet

### Strep-Tag II antibody (orb499988)

<b>Description</b>	Strep-Tag II antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	All
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	KLH conjugated synthetic peptide Strep-Tag
<b>Target</b>	Strep-Tag II
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	WB=1:1000-5000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot: Lane 1: Strep-Tagged Fusion ...