

Product Datasheet

MLKL (Phospho-Ser358) antibody (orb499999)

Description	MLKL (Phospho-Ser358) antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised monomethylpeptide derived from human MLKL around the phosphorylation site of Ser358 (QT(p-S)MS)
Target	phospho-MLKL (Ser358)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	54 kDa
Uniprot ID	Q9D2Y4
Expiration Date	12 months from date of receipt.