

Product Datasheet

Rabbit IgG(Fc Specific,Heavy Chain) antibody (orb500025)

Description	Rabbit IgG(Fc Specific,Heavy Chain) antibody
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	rabbit IgG
Target	Rabbit IgG(Fc specific,Heavy Chain)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	55 kDa
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.