


Product Datasheet

Human BMP-4 ELISA Kit (orb50010)

Description	Human BMP-4 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S293-R408
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum is 0-110 pg/ml
Sample Types	cell culture supernatants, cell lysates and serum.
MW	46555 MW
Uniprot ID	P12644
Sensitivity	4 pg/ml
Expiration Date	Please enquire.



Human BMP-4 PicoKine ELISA Kit standard ...