


Product Datasheet

Mouse Eotaxin-1 ELISA Kit (orb50016)

Description	Mouse Eotaxin-1 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: H24-P97
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is around 1 ng/ml
Sample Types	cell culture supernatants, serum and plasma (citrate, EDTA).
MW	10893 MW
Uniprot ID	P48298
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Mouse
Eotaxin
PicoKine
ELISA Kit
standar...