


Product Datasheet

Human IFN Gamma/IFNG/Interferon Gamma ELISA Kit (orb50030)

Description	Human IFN Gamma/IFNG/Interferon Gamma ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q24-Q166
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 15.6 pg/ml
Sample Types	cell culture supernatants and serum.
MW	19348 MW
Uniprot ID	P01579
Sensitivity	2 pg/ml
Expiration Date	Please enquire.



Human IFN gamma PicoKine ELISA Kit stand...