


Product Datasheet

Mouse IFN Gamma / IFNG / Interferon gamma ELISA Kit (orb50031)

Description	Mouse IFN Gamma / IFNG / Interferon gamma ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: H23-C155
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 31.2 pg/ml
Sample Types	cell culture supernatants and serum.
MW	17907 MW
Uniprot ID	P01580
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Mouse IFN gamma PicoKine ELISA Kit stand...