


Product Datasheet

Nano-Tag (9) antibody (orb500329)

Description	Nano-Tag (9) antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide Nano-Tag (9)
Target	Nano-Tag (9)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	WB=1:1000-10000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Nano-Tagged Fusion P...