


Product Datasheet

Rat IGF-1 ELISA Kit (orb50033)

Description	Rat IGF-1 ELISA Kit
Reactivity	Rat
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G49-A118
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:500, concentration in serum and plasma is around 142-340 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	17079 MW
Uniprot ID	P08025
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Rat IGF-1
PicoKine
ELISA Kit
standard
cu...