


Product Datasheet

Human IGFBP-1 ELISA Kit (orb50036)

Description	Human IGFBP-1 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A26-N259
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:20, concentration in serum and plasma is 159-75205 pg/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	27904 MW
Uniprot ID	P08833
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Human
IGFBP-1
PicoKine
ELISA Kit
standar...