


Product Datasheet

Maltose Binding Protein antibody (orb500360)

Description	Maltose Binding Protein antibody
Species/Host	Mouse
Reactivity	E. coli, Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant MBP Protein of Escherichia coli (strain K12)
Target	MBP tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	44 kDa
Uniprot ID	POAEX9
Dilution Range	WB=1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: MBP tagged cytochrom...