


## Product Datasheet

### Mouse IGFBP-1 ELISA Kit (orb50037)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Mouse IGFBP-1 ELISA Kit   |
| <b>Reactivity</b>          | Mouse   |
| <b>Range</b>               | 31.2 pg/ml - 2,000 pg/ml  |
| <b>Tested Applications</b> | ELISA   |
| <b>Immunogen</b>           | Expression system for standard: NS0; Immunogen sequence: A26-N272   |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Dilution ratio of 1:10, concentration in serum and plasma is 2737-4530 pg/ml  |
| <b>Sample Types</b>        | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).  |
| <b>MW</b>                  | 29570 MW  |
| <b>Uniprot ID</b>          | <b>P47876</b>   |
| <b>Sensitivity</b>         | 3 pg/ml   |
| <b>Expiration Date</b>     | Please enquire.   |



Mouse  
IGFBP-1  
PicoKine  
ELISA Kit  
standar...