


Product Datasheet

Human IL-1 Beta/IL-1F2/IL1B ELISA Kit (orb50042)

Description	Human IL-1 Beta/IL-1F2/IL1B ELISA Kit
Reactivity	Human
Range	3.9 pg/ml - 250 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A117-S269
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 12.5 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	30748 MW
Uniprot ID	P01584
Sensitivity	0.15 pg/ml
Expiration Date	Please enquire.



Human IL-1
beta
PicoKine
ELISA Kit
stand...