



## Product Datasheet

### GST Tag antibody (orb500513)

<b>Description</b>	GST Tag antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	All
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant GST Protein
<b>Target</b>	GST tag
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	WB=1:30000-100000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot:GST-Tagged fusion protein O...



Western blot:Recombinant GST Protein;Pri...