


Product Datasheet

Human IL-6/Interleukin-6 ELISA Kit (orb50052)

Description	Human IL-6/Interleukin-6 ELISA Kit
Reactivity	Human
Range	4.69 pg/ml - 300 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: P29-M212
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 4.69 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	23,718 MW
Uniprot ID	P05231
Sensitivity	0.3 pg/ml
Expiration Date	Please enquire.



Human IL-6
PicoKine
ELISA Kit
standard
C...