

Product Datasheet

Beta 2 Microglobulin antibody (orb500534)

Description	Beta 2 Microglobulin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human Beta-2-Microglobulin protein
Target	beta 2 Microglobulin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13/27 kDa
Uniprot ID	P61769
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.