


Product Datasheet

Human CCL22/MDC ELISA Kit (orb50072)

Description	Human CCL22/MDC ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G25-Q93
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 500 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10625 MW
Uniprot ID	000626
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Human
CCL22/MDC
PicoKine
ELISA Kit
stand...