


Product Datasheet

Human MIP-1Alpha/CCL3 ELISA Kit (orb50074)

Description	Human MIP-1Alpha/CCL3 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A27-A92
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	10085 MW
Uniprot ID	AONOR1
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human MIP-1 alpha PicoKine ELISA Kit sta...