


Product Datasheet

Human MIP-3 Alpha / CCL20 / MIP3a ELISA Kit (orb50075)

Description	Human MIP-3 Alpha / CCL20 / MIP3a ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A27-M96
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 0-171 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10762 MW
Uniprot ID	P78556
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Human MIP-3
alpha/CCL20
PicoKine
ELISA K...