


Product Datasheet

PI 3 Kinase P85 Alpha antibody (orb500765)


Description	PI 3 Kinase P85 Alpha antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombinant mouse PI3K p85 Protein (1-110/724 aa)
Target	PIK3R1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	80 kDa
Uniprot ID	P27986
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



ICC analysis of SH-SY5Y cell using PI 3 ...



IHC-P analysis of rat skeletal muscle us...



IHC-P analysis of Mouse brain using PI3K...