



## Product Datasheet

### Beta Tubulin antibody (orb500779)


|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Beta Tubulin antibody                                                                                                                   |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | FC, IHC-P, WB                                                                                                                           |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human Beta tubulin (61-160/444 aa)                                                        |
| <b>Target</b>              | Beta tubulin                                                                                                                            |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                                                                   |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                  |
| <b>Concentration</b>       | 1mg/ml                                                                                                                                  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                                              |
| <b>MW</b>                  | 55 kDa                                                                                                                                  |
| <b>Uniprot ID</b>          | <b>P07437</b>                                                                                                                           |
| <b>Dilution Range</b>      | WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1µg/Test, IHC-F=1:100-500, IF=1:100-500                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |



Flow cytometry of HL-60.Primary Antibody...



IHC-P analysis of mouse cerebellum using...



Western blot:Cerebellum (Mouse);Cerebell...