



Product Datasheet

Histone H3 antibody (orb500825)


Description	Histone H3 antibody
Species/Host	Mouse
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Histone H3
Target	Histone H3 (Nuclear Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	15 kDa
CAS Number	Histone H3.1
Uniprot ID	P68431
Dilution Range	WB=1:20000, IHC-P=1:100-500, IF=1:100-500, IHC-F=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of mouse liver using Hist...



IHC-P analysis of rat testis using Histo...



IHC-P analysis of human gastric carcinom...