


## Product Datasheet

### Human PDGF-AB ELISA Kit (orb50095)

<b>Description</b>	Human PDGF-AB ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.2 pg/ml - 2,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: E.coli; Immunogen sequence: (S87-T211) + (S82-T190)
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	For serum, dilution ratio is 1:20 and concentration is 5432-34260 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 34-732 pg/ml
<b>Sample Types</b>	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
<b>MW</b>	27283 MW
<b>Uniprot ID</b>	<b>P04085, P01127</b>
<b>Sensitivity</b>	2 pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
PDGF-AB  
PicoKine  
ELISA Kit  
standar...