

Product Datasheet

pCRMP2 antibody (Biotin) (orb501048)

Description	pCRMP2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human CRMP2 around the phosphorylation site of Thr514+Ser518
Target	Phospho-CRMP2 (Thr514+Ser518)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16555
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500
Expiration Date	12 months from date of receipt.