

Product Datasheet

RACGAP1(phospho-Ser387) antibody (Biotin) (orb501057)

Description	RACGAP1(phospho-Ser387) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human RACGAP1 around the phosphorylation site of Ser387
Target	phospho-RACGAP1 (Ser387)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H0H5
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.