

Product Datasheet

NBS1 p95 (phospho-Ser343) antibody (Biotin) (orb501068)

Description	Rabbit polyclonal antibody against NBS1 p95 (phospho-Ser343) conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NBS1 around the phosphorylation site of Ser343 SL(p-S)QG
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85kDa
Uniprot ID	O60934
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.