


Product Datasheet

Human TIE2/TEK ELISA Kit (orb50107)

Description	Human TIE2/TEK ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A23-K745
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is 14-90 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	125830 MW
Uniprot ID	Q02763
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Human TIE2
PicoKine
ELISA Kit
standard c...