

Product Datasheet

SCNN1B (phospho-Ser633) antibody (Biotin) (orb501079)

Description	SCNN1B (phospho-Ser633) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SCNN1B around the phosphorylation site of Ser633
Target	phospho-SCNN1B (Ser633)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51168
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.