


Product Datasheet

Human TIMP-3 ELISA Kit (orb50110)

Description	Human TIMP-3 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C24-P211
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 600-1200 pg/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and saliva.
MW	24145 MW
Uniprot ID	P35625
Sensitivity	2 pg/ml
Expiration Date	Please enquire.



Human TIMP-3 PicoKine ELISA Kit standard...