

Product Datasheet

PRKD1 (phospho-Ser738) antibody (Biotin) (orb501125)

Description	Rabbit polyclonal antibody against PRKD1 (phospho-Ser738) conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PRKD1 around the phosphorylation site of Ser738 EK(p-S)FR
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	102kDa
Uniprot ID	Q15139
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.