


Product Datasheet

Human sTNFsR II ELISA Kit (orb50114)

Description	Human sTNFsR II ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L23-D257
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:20, concentration in serum and plasma is 802-3804 pg/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
MW	48291 MW
Uniprot ID	P20333
Sensitivity	2 pg/ml
Expiration Date	Please enquire.



Human
sTNFsR II
PicoKine
ELISA Kit
stand...