


## Product Datasheet

### Human VEGFR1/sFLT1 ELISA Kit (orb50122)

<b>Description</b>	Human VEGFR1/sFLT1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	156 pg/ml - 10,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: NS0; Immunogen sequence: S27-H687
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	For serum, dilution ratio is 1:1 and concentration is 3.6-5 ng/ml. For plasma, dilution ratio is 1:1 and concentration is less than the lowest standard, 156 pg/ml
<b>Sample Types</b>	cell culture supernatants, serum and plasma (EDTA).
<b>MW</b>	150769 MW
<b>Uniprot ID</b>	<b>P17948</b>
<b>Sensitivity</b>	30 pg/ml
<b>Expiration Date</b>	Please enquire.



Human VEGFR1/sFLT1 PicoKine ELISA Kit st...