


Product Datasheet

Human E-Cadherin/CDH1/Cadherin-1 ELISA Kit (orb50124)

Description	Human E-Cadherin/CDH1/Cadherin-1 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D155-I707
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is 17-71 ng/ml
Sample Types	cell culture supernatants, serum, plasma (heparin) , saliva and urine.
MW	97456 MW
Uniprot ID	P12830
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Human E-Cadherin PicoKine ELISA Kit stan...