


Product Datasheet

Human Resistin ELISA Kit (orb50128)

Description	Human Resistin ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: K19-P108
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is around 20 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11419 MW
Uniprot ID	Q9HD89
Sensitivity	3 pg/ml
Expiration Date	Please enquire.



Human Resistin ELISA Kit standard curve