

Product Datasheet

GSK3 beta (phospho-Tyr216) antibody (Biotin) (orb501287)

| | |
|----------------------------|---|
| Description | GSK3 beta (phospho-Tyr216) antibody (Biotin) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Biotin |
| Tested Applications | IHC-P, WB |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human GSK-3 Beta around the phosphorylation site of Tyr216 |
| Target | phospho-GSK3 Beta (Tyr216) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P49841 |
| Dilution Range | WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |