


Product Datasheet

Human VEGF-C ELISA Kit (orb50131)

Description	Human VEGF-C ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T103-R227
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 500-1250 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	46883 MW
Uniprot ID	P49767
Sensitivity	3 pg/ml
Expiration Date	Please enquire.



Human VEGF-C PicoKine ELISA Kit standard...