


Product Datasheet

Rat IL-18/IL-1F4/Interleukin-18 ELISA Kit (orb50132)

Description	Rat IL-18/IL-1F4/Interleukin-18 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: H37-S194
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 15.6 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	22303 MW
Uniprot ID	P97636
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Rat IL-18
PicoKine
ELISA Kit
standard
cu...