
Product Datasheet

**Rb (phospho-Thr821) antibody (Biotin)
(orb501381)**

Description

Rabbit polyclonal antibody against Rb

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Retinoblastoma associated protein around the phosphorylation site of Thr821
Target	RB1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	102kDa
Uniprot ID	P06400
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.