

Product Datasheet

NDRG1 (phospho-Thr346) antibody (Biotin) (orb501422)

Description	NDRG1 (phospho-Thr346) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NDRG1 around the phosphorylation site of Thr346
Target	Phospho-NDRG1 (Thr346)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q92597
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.