

## Product Datasheet

# NFKB1 (phospho-Ser927) antibody (Biotin) (orb501424)

<b>Description</b>	NFKB1 (phospho-Ser927) antibody (Biotin)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Rat
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, IHC-Fr, IHC-P
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from rat NFKB1 around the phosphorylation site of Ser927
<b>Target</b>	Phospho-NFKB1 (Ser927)
<b>Preservatives</b>	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P19838</b>
<b>Dilution Range</b>	IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
<b>Expiration Date</b>	12 months from date of receipt.