

Product Datasheet

LKB1 (phospho-Thr189) antibody (Biotin) (orb501455)

Description	LKB1 (phospho-Thr189) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human LKB1 around the phosphorylation site of Thr189
Target	Phospho-LKB1 (Thr189)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15831
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.