


Product Datasheet

Human MIF ELISA Kit (orb50147)

Description	Human MIF ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: P2-A115
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is 5-13 ng/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and human milk.
MW	12476 MW
Uniprot ID	I4AY87
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human MIF ELISA Kit standard curve