


Product Datasheet

Human IL-18/IL-1F4/Interleukin-18 ELISA Kit (orb50153)

Description	Human IL-18/IL-1F4/Interleukin-18 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Y33-D197
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 67-429 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	22326 MW
Uniprot ID	Q14116
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Human IL-18 PicoKine ELISA Kit standard ...