

---

## Product Datasheet

**AMPK Beta 1 + AMPK Beta 2 (phospho-Ser182 antibody (Biotin) (orb501554)**

**Description**

AMPK Beta 1 + AMPK Beta 2

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Mouse, Rat
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, IHC-Fr, IHC-P, WB
<b>Immunogen</b>	KLH conjugated synthesised phosphopeptide derived from human AMPK beta 1 around the phosphorylation site of Ser182
<b>Target</b>	phospho-AMPK beta 1 + AMPK beta 2 (Ser182/Ser184)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q9Y478</a>
<b>Dilution Range</b>	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
<b>Expiration Date</b>	12 months from date of receipt.