

Product Datasheet

ATP citrate lyase (phospho-Ser455) antibody (Biotin) (orb501562)

Description	ATP citrate lyase (phospho-Ser455) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ACLY around the phosphorylation site of Ser455
Target	phospho-ACLY (Ser455)
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P53396
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.