

---

## Product Datasheet

**MAPK14 (phospho-Thr180/Tyr182) antibody  
(Biotin) (orb501616)**

**Description**

MAPK14 (phospho-Thr180/Tyr182)

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	IHC-Fr, IHC-P, WB
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human p38 MAPK around the phosphorylation site of Thr180/Tyr182
<b>Target</b>	Phospho-P38 MAPK (Thr180 + Tyr182)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q16539</a>
<b>Dilution Range</b>	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500
<b>Expiration Date</b>	12 months from date of receipt.